

# Add collection event0

After Arms are added to Epochs, next step is to add *Collection events* to an arm. To add a collection event, select the desired Arm from the tree and select *Add collection event* button.

Study details

- 10600
  - Pre-Registration
    - Pre-Registration
  - Remission Induction Therapy
  - Midostaurin

EDIT ARM

Label

Pre-Registration

Event label	Clinical status	Calendar event point
-------------	-----------------	----------------------

Submit Cancel Add collection event

Save study Synchronize study

Add collection event screen appears. Enter an *Event label*, from the *Clinical status* drop-down select a value. Enter the *Calendar event point*, in *Clinical diagnosis* filed, start typing the clinical diagnosis and system will suggest the list of values that match the entered text. From the suggestions, select a clinical diagnosis and select *Submit*.

Study details

- 10600
  - Pre-Registration
    - Pre-Registration
  - Remission Induction Therapy
  - Midostaurin

ADD COLLECTION EVENT

Event label \* : Screening

Clinical status \* : Not Specified

Calendar event point \* : 2.0 Days

Clinical diagnosis \* : [M]Medulloblastoma NOS

Submit Cancel Add Specimen Requirement

Save study Synchronize study

The collection event is added to the selected arm and system takes you to a page where collection events for an arm are displayed. To add another collection event, select *Add collection event* button in this screen and follow the steps that are followed above. Any number of collection events may be added to an arm.

Study details

- 10600
  - Pre-Registration
    - Pre-Registration**
    - Screening
  - Remission Induction Therapy
  - Midostaurin

EDIT ARM

Label

Pre-Registration

Event label	Clinical status	Calendar event point
Screening	Not Specified	2.0

Submit Cancel Add collection event

Save study Synchronize study