

Add alternate specimen0

When CRA collects specimens from a patient, in absence of a desired specimen, he/she may collect any other specimen that is added as alternate to it. To add an alternate specimen to a specimen requirement, select the desired specimen requirement from the tree structure and navigate to Alternate specimens section and select *Add* button.

Study details

- 10600
 - Pre-Registration
 - Pre-Registration
 - Screening
 - 5 ml Bone marrow Aspirate (Green Top)**
 - Remission Induction Therapy
 - Midostaurin

ADD SPECIMEN REQUIREMENT

Title * 5 ml Bone marrow Aspirate (Green Top)

Class * Fluid Type * Not Specified

Anatomical site * Not Specified Laterality of site * Not Specified

Quantity * 5.0 ml Pathological status * Not Specified

Tube type * Not Specified Concentration

Preparation instructions Draw bone marrow aspirate in a green top tube (heparinized). Label the tube with the patients CALGB patient ID number,

Alternate specimen(s) **Add** Remove

SHIPPING REQUIREMENTS

Shipment site * WASHU Repository Shipping method * Cold Pack

Shipping Instructions Samples must be shipped on the same day they are obtained. Send the samples at room temperature, any day of the week, via overnight

Submit **Cancel** **Add equivalent specimens**

System displays Add alternate specimen screen for the selected specimen requirement. Add the details related to alternate specimen and select *Submit*.

Study details

- 10600
 - Pre-Registration
 - Pre-Registration
 - Screening
 - 5 ml Bone marrow Aspirate (Green Top)
 - Remission Induction Therapy
 - Midostaurin

ADD ALTERNATE SPECIMEN

Title *

Class * **Type ***

Anatomical site * **Laterality of site ***

Quantity * **Pathological status ***

Tube type * **Concentration**

Preparation instructions

Alternate specimen(s)

SHIPPING REQUIREMENTS

Shipment site * **Shipping method ***

Shipping Instructions

System adds an alternate specimen to the selected specimen requirement and tree structure is updated to display the alternate specimen

Specimen requirement added successfully

Study details

- 10600
 - Pre-Registration
 - Pre-Registration
 - Screening
 - 5 ml Bone marrow Aspirate (Green Top)
 - Bone Marrow Biopsy Section
- Remission Induction Therapy
- Midostaurin

ADD ALTERNATE SPECIMEN

Title *

Class * Type *

Anatomical site * Laterality of site *

Quantity * gm Pathological status *

Tube type * Concentration

Preparation instructions

Alternate specimen(s)

SHIPPING REQUIREMENTS

Shipment site * Shipping method *

Shipping Instructions

After selecting submit, select *Save study* button. You may add any number of alternate specimens to a specimen requirement. All of them are displayed in the tree structure as well as in the Alternate specimens section of that specimen requirement.

Study details

- 10600
 - Pre-Registration
 - Pre-Registration
 - Screening
 - 5 ml Bone marrow Aspirate (Green Top)**
 - Bone Marrow Biopsy Section
 - Remission Induction Therapy
 - Midostaurin

Save study Synchronize study

ADD SPECIMEN REQUIREMENT

Title *

Class *

Anatomical site *

Quantity * ml

Tube type *

Type *

Laterality of site *

Pathological status *

Concentration

Preparation instructions

Alternate specimen(s)

SHIPPING REQUIREMENTS

Shipment site * Shipping method *

Shipping Instructions

Only when you select the *Class* as Molecular, system activates the concentration field. Enter a value for concentration and follow the other steps mentioned above to complete adding a specimen requirement