

Add epoch0

To add an epoch to a study, select *Add epoch* button

The screenshot shows the 'Study details' form for study 10600. The form includes fields for Alliance study number (10600), Title (10600 - Treatment study), Start date (07/01/2012), Author (Admin BMS), and Clinical diagnosis (Abnormality of neurotensine secretion (disorde [M]Dermatofibroma NOS). The 'Add epoch' button is highlighted with an orange box.

Add epoch screen appears. It is initialized with a single row.

The screenshot shows the 'ADD EPOCH(S)' screen. It features a table with two columns: 'Label' and 'Sequence'. The table contains one row with the label 'Pre-Registration' and a sequence value of 1. The 'Submit' button is highlighted with an orange box.

<input type="checkbox"/>	Label	Sequence
<input type="checkbox"/>	Pre-Registration	1

Enter Label and Sequence for the epoch and select *Submit*. To add more epochs, select *Add more* button. A new row is added every time you select the *Add more* button.

Study details

- 10603
 - Pre-Registration

ADD EPOCH(S) Add more Delete

	Label	Sequence
<input type="checkbox"/>	Pre-Registration	1
<input type="checkbox"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/>	<input type="text"/>	<input type="text"/>

Submit Cancel

Save study
Synchronize study

Enter the *Label* and *Sequence* for the new rows and select *Submit* button. Note that for a study label of an epoch, sequence are unique. Hence, no two epochs may have the same *Label* or *Sequence*.

Study details

- 10603
 - Pre-Registration
 - Remission Induction Therapy
 - Midostaurin

ADD EPOCH(S) Add more Delete

	Label	Sequence
<input type="checkbox"/>	Pre-Registration	1
<input type="checkbox"/>	Remission Induction Therapy	2
<input type="checkbox"/>	Midostaurin	3

Submit Cancel

Save study
Synchronize study

Once you select *Submit*, labels of the epochs are displayed in the tree structure.